Ex.A.O.

Capteur Foxy/Tooxy/AirNeXT

Data logging

Sensor

Réf :
482 056

Français – p 1

English – p 5

Version : 8110

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0 825 563 563*

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1 General Information

This sensor measures conductivity in ExAO.
It has a conductivity meter probe CTA and offers automatic temperature off-setting.

2 Utilisation

The sensor operates with the Foxy®/AirNeXT® console and can be used:

- In standalone mode, displaying the results on the sensor display,
- In computer connected mode, displaying the results on the computer.

2.1 Connection to the Foxy®/AirNeXT® console

The sensor is connected using the connector located on the unit:

- Position the sensor so it is aligned with the selected measuring channel,
- Engage the sensor on the guides in the console,
- Insert firmly to the end.

2.2 Type of measurement

The sensor, fitted with the probe, measures the conductivity of a solution.

2.3 Operation in standalone mode

Starting-up:


- Connect the sensor to the console,
- Power on the console.

- **Calibration:**

First, press and briefly hold down the zero-setting button (fuchsia button); the arrow above "Solution" will flash. ▼

Immerse the probe in a beaker containing a reference standard solution, rotate the potentiometer to display the reference standard value then briefly press the same fuchsia adjustment button; the arrow above "Solution" will go out.

The conductivity meter probe is now calibrated.

Display:

A conductivity value is displayed on the conductivity meter sensor screen.

Type of measurement/Calibration:

Changes of calibration are made using the green button provided for this purpose.


2.4 Operation in connected mode

Starting-up:


- Connect the sensor to the console,
- Power on the console,
- Connect the console to the computer.

Calibration:

An “Adjustment” tab on the acquisition software is used to calibrate the sensor.

Click on the “Start” button after leaving the probe in a beaker containing a reference standard solution and wait for the value to stabilise.
A red square flashes around the “Adjustment” button when the value is stabilising.
When the square turns green the value has stabilised. Then click on “Adjust” after entering the reference value to confirm the calibration.

Display:

The conductivity meter sensor display is inactive, measurement are displayed directly on the computer screen via the acquisition software.

Type of measurement/Calibration:

Changes of calibration are made using the acquisition software and the sensor button is deactivated.

3 Technical Features

Calibrations	2mS.cm ⁻¹	20mS.cm ⁻¹
Resolution	2μS.cm ⁻¹	20μS.cm ⁻¹
Precision	1%	1%
Temperature off-set	Automatic	0 to 50°C

4 After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.